

DENTAL CARE PRODUCTS

SUPER 2000 MASSAGE TOOTH BRUSH

It is estimated that 75% of adults have gum disease and 80% of kids over the age of 8 have some form of it. Gum disease is the second most preventable disease next to the common cold. The best solution to prevent gum disease is proper brushing.

Conventional toothbrushes require skill and longer time to achieve proper brushing. The *Super 2000* Massage Tooth brush is designed to assist adults and kids to brush easier, better, quicker and more effectively.

With the high speed vibration of 2000 stokes per minutes, *Super 2000* toothbrush helps you to brush front, back and chewing surfaces effortlessly. Simply move the device gently across all areas and your teeth will be clean and fresh after each application. *Super 2000* also massages your front and back gums to relax them and prevent disease.

Super 2000 Massage Toothbrush is operated by one "AA" battery with easy installation. It is safe and user friendly. The package consists of 5 brushes and 2 batteries for the entire family uses. It is the most complete and economical electric toothbrush available on the market today.

Super 2000 Massage Toothbrush cleans teeth, removes plaque, prevents gum diseases and relieves tension by massaging your entire mouth. No wonder many dentists around the world recommend their patients to use it daily.

BATTERY INSTALLATION

1. Turn leftward to open the back cover.
2. Install one "AA" battery properly.
3. Turn rightward and close the back cover.

DIRECTIONS FOR USE

1. Put on brush head.
2. Turn the bottom portion to the "on" position for tooth brushing.
3. Move brush gently across the entire surface of teeth and along gumline to clean.
4. For better results, brush teeth twice daily at bedtime and in the morning.
5. After you finish brushing, turn the unit back to "off" position.